# Peri-Strips Dry®

Staple Line Reinforcement

Klammernaht-Verstärkung

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Steam heat sterilized

NaOH

Product treated with sodium hydroxide.

BOVINE

Product derived from USDA-inspected cattle.

Rx Only

CAUTION: Federal (USA) law restricts this device to sale by, or on the order of, a physician.

REF

Catalog number

Use by

LOT

Lot number

Compatible stapler models

## DESCRIPTION:

Peri-Strips Dry® Staple Line Reinforcement (Peri-Strips Dry) is a biologic tissue prepared from bovine pericardium procured from cattle originating in the United States. Peri-Strips Dry has been chemically treated using glutaraldehyde, ethanol, and propylene oxide. Peri-Strips Dry is treated with 1 molar sodium hydroxide for 60 - 75 minutes at 20 - 25°C, a method recommended by international experts for reducing or inactivating prions, the agent suspected of causing Transmissible Spongiform Encephalopathies (TSE)/Bovine Spongiform Encephalopathies (BSE).

The use of Peri-Strips Dry requires two primary components: The Peri-Strips Dry tissue assembly (plastic sheath, foam, and two Peri-Strips Dry strips and PSD Gel. One (1) tube of PSD Gel (Gel) is provided for every two (2) pouches of Peri-Strips Dry. The Gel is used to create a temporary bond between the Peri-Strips Dry strips and the surgical stapler forks until the stapler is positioned and fired. Each Peri-Strips Dry tissue assembly and each Gel tube is packaged sterile in a separate pouch.

## INDICATIONS:

To reinforce staple lines during lung and bronchus resections including: pneumonectomy, pneumoreduction, pneumectomy, segmentectomies (segmented resections), wedge resections, blebectomies, lobectomies, bullectomies, bronchial resections and other lung incisions and excisions of lung and bronchus. This product may also be used for reinforcement of the gastric staple line during the bariatric surgical procedures of gastric banding and gastric bypass in the United States, European Union, Canada, and Australia.

Peri-Strips Dry® is intended to be used with the surgical stapler models indicated on the package labeling. Use with any other stapler models is not recommended.

## CONTRAINDICATIONS:

The use of Peri-Strips Dry is contraindicated in patients with known sensitivity to bovine material.

## ADVERSE REACTIONS:

As with any surgical procedure, adverse reactions are possible and include but are not limited to: infection, rejection, erosion, and allergic reaction.

## WARNINGS:

Do not re-sterilize.

Do not use product if there is damage to the pouch or seals.

Ensure the stapler forks are completely covered with the strips, or the strips will not be below the staples after firing and may result in inadequate repair.

Peri-Strips Dry is not designed, sold, or intended for use except as indicated; doing so may result in surgical complications.

Synovis products differ; substitution of one product for another product may be harmful to the patient.

Glutaraldehyde-treated bovine pericardium may undergo accelerated calcific infiltration in patients with high calcium metabolic activity (e.g., children). This may not be a concern where the patch is exposed to systolic pressures.

When used to correct simple complete transposition of the great arteries with pericardium augmentation of the pulmonary venous channel, bovine pericardium has been reported to demonstrate calcification, inflammation, and formation of fibrous tissue which obstructed pulmonary venous flow.

#### CAUTIONS:

Do not apply excess pressure to the tissue assembly when removing the stapler from the tissue assembly, or the strips may not adhere to the stapler properly.

Do not get the strips wet before applying Gel, or the strips may not adhere to the stapler properly.

Ensure the anvil and cartridge sides of the tissue assembly are on the corresponding stapler forks or the strips may not adhere to the stapler properly.

Final tissue compression, including Peri-Strips Dry, must meet the range specified by the stapler manufacturer; this is especially important if staple firings are overlapped. Peri-Strips Dry increases the total thickness of the area stapled by 0.5 mm - 1.5 mm (0.02" - 0.06").

Follow Instructions for Use supplied by the stapler manufacturer. Do not use Peri-Strips Dry contrary to the stapler manufacturer's instructions.

## INSTRUCTIONS FOR USE:

### I. LOADING STAPLER

**Note:** Each model of Peri-Strips Dry has been designed specifically for the stapler models indicated on the label; verify that the correct model of Peri-Strips Dry has been selected.

**Note:** Loading technique for Peri-Strips Dry varies. Please follow appropriate technique as indicated below.

**Instruction for Peri-Strips Dry with the following codes: PSD 4502/6-EN, PSD 4502/6-ETSN, PSD 6002/6 - UN (when used with non-cutting stapler), PSD 6002/6-ECH, PSD 6002/6-U-DST, PSD 8002/6-U-DST, PSD 9002/6- PIN, PSD 9002/6-UN, PSD 8002/6-UN, PSD 7502/6-EN, PSD-5502/6 EN, PSD 10002/6- EN, PSD 7502/6-EN and PSD 5502-EN (See figure 1)**


FIGURE 1

- 1a. Peel open the outer Peri-Strips Dry pouch, aseptically remove the inner pouch; the inner pouch may be placed in the sterile field. Inspect the pouch. Do not use if the pouch is damaged or if the seals are not intact.
- 2a. Open the inner pouch. Use sterile atraumatic techniques to remove the Peri-Strips Dry tissue assembly .
- 3a. Inspect the Gel pouch. Do not use if the pouch is damaged or if the seals are not intact. Peel open the Gel pouch, aseptically remove the Gel and place in the sterile field.

- 4a. Break off the tip of the Gel tube (See figure 2).


FIGURE 2

- 5a. Attach plastic tip to the Gel tube by pushing it down over the top of the Gel tube (See figure 3).

FIGURE 3

- 6a. Apply a continuous thin Gel bead onto each strip (See figure 4).

**Caution:** Do not get the strips wet before applying Gel or the strips may not adhere to the stapler properly.


FIGURE 4

- 7a. Ensure the forks of the stapler are clean and dry before inserting the Peri-Strips Dry tissue assembly.

**Caution:** Final tissue compression, including Peri-Strips Dry, must meet the range specified by the stapler manufacturer; this is especially important if staple firings are overlapped. Peri-Strips Dry increases the total thickness of the area stapled by 0.5 mm - 1.5 mm (0.02" - 0.06").

- 8a. Identify the anvil (ANV) and cartridge (CART) sides of the Peri-Strips Dry tissue assembly. Position the open stapler onto the tissue assembly (See figure 5).

**Caution:** Ensure the anvil and cartridge sides of the tissue assembly are on the corresponding stapler forks or the strips may not adhere to the stapler properly.


FIGURE 5

- 9a. Close the stapler.

- 10a. Remove plastic sheath, leaving the foam spacer between the stapler forks. Discard the plastic sheath and remove any excess gel.

**Note:** For optimal adherence of Peri-Strips Dry onto the stapler jaws, wait at least 15 seconds before opening the jaws and removing the foam spacer.

**Note:** For optimal results, the prepared Peri-Strip Dry/stapler should be used within 60 minutes.


FIGURE 6

11 a. Just prior to using the stapler, open stapler and remove foam spacers (See figure 7). Visually inspect strips to ensure each Peri-Strips Dry strip lies flat and fully covers the staple lines.


FIGURE 7

Instructions for Peri-Strips Dry with the following product codes: PSD 3002/6-UN, PSD 4502/6-UN and PSD 6002/6-UN (when used with cutting stapler) (See figure 8)


FIGURE 8

- 1b. Peel open the outer Peri-Strips Dry pouch, aseptically remove the inner pouch; the inner pouch may be placed in the sterile field. Inspect the pouch. Do not use if the pouch is damaged or if the seals are not intact.
- 2b. Open the inner pouch. Use sterile atraumatic techniques to remove the Peri-Strips Dry tissue assembly.
- 3b. Inspect the Gel pouch. Do not use if the pouch is damaged or if the seals are not intact. Peel open the Gel pouch, aseptically remove the Gel and place in the sterile field.
- 4b. Break off the tip of the Gel tube (See figure 9).


FIGURE 9

- 5b. Attach plastic tip to the Gel tube by pushing it down over the top of the Gel tube (See figure 10).

FIGURE 10

- 6b. Apply a continuous thin Gel bead onto each strip (See figure 11).

**Caution:** Do not get the strips wet before applying Gel or the strips may not adhere to the stapler properly.


FIGURE 11

- 7b. Ensure the forks of the stapler loading unit are clean and dry before inserting the Peri-Strips Dry tissue assembly.

**Caution:** Final tissue compression, including Peri-Strips Dry, must meet the range specified by the stapler manufacturer; this is especially important if staple firings are overlapped. Peri-Strips Dry increases the total thickness of the area stapled by 0.5 mm - 1.5 mm (0.02" - 0.06").

- 8b. Identify the anvil (ANV) and cartridge (CART) sides of the Peri-Strips Dry tissue assembly. Holding the cartridge side up, slide stapler loading unit onto the tissue assembly, being careful to ensure that the end of the yellow guard is positioned above the assembly clip (See figure 12).


FIGURE 12

**Caution:** Ensure the anvil and cartridge sides of the tissue assembly are on the corresponding fork of the stapler loading unit or the strips may not adhere to the stapler properly.

- 9b. While applying light pressure to both sides of the stapler forks, slide assembly clip over the stapler loading unit, ensuring that both sides of the stapler are held firmly in place by the assembly clip and that, on the cartridge side, the yellow guard is positioned between the assembly clip (See figure 13).

**Note:** Peri-Strips Dry can be used immediately or allowed to remain between the forks of the stapler loading unit until surgery.

**Note:** For optimal results, the prepared Peri-Strips Dry/stapler should be used within 60 minutes.


FIGURE 13

- 10b. Immediately before using, slide assembly clip off stapler forks. Remove stapler loading unit (See figure 14).
- 11b. Visually inspect strips to ensure each Peri-Strips Dry strip lies flat and fully covers the staple lines.


FIGURE 14

## II. IMPLANTING PERI-STRIPS DRY

12. Follow *Instructions For Use* supplied by the stapler manufacturer.
13. If necessary, cut the end of the Peri-Strips Dry to remove the dissected tissue (See figure 15).


FIGURE 15

14. Discard any opened Peri-Strips Dry and Gel pouches. These cannot be re-sterilized or re-used.

## DISCLAIMER OF WARRANTIES

Synovis Surgical Innovations (SSI), a division of Synovis Life Technologies, Inc., warrants that reasonable care has been used in the manufacture of this device. This warranty is exclusive and in lieu of all other warranties whether expressed, implied, written or oral, including, but not limited to, any implied warranties of merchantability or fitness. As a result of biological differences in individuals, no product is 100% effective under all circumstances. Because of this fact and since SSI has no control over the conditions under which the device is used, diagnosis of the patient, methods of administration or its handling after it leaves its possession, SSI does not warrant either a good effect or against an ill effect following its use. The manufacturer shall not be liable for any incidental or consequential loss, damage or expense arising directly or indirectly from the use of this device. SSI will replace any device which is defective at the time of shipment. No representative of SSI may change any of the foregoing or assume any additional liability or responsibility in connection with this device.

## REFERENCES:

World Health Organization. Public health issues related to animal and human spongiform encephalopathies: memorandum from a WHO meeting. *Bul WHO*. 1992; 70(2); 183 - 190.

EEC Document. Guidelines for minimizing the risk of transmitting agents causing spongiform encephalopathy via medicinal products. *Biologicals*. 1992; 20; 155 - 158.

DiMartino A. et al. Inactivation of the scrapie agent in a scaled-down procedure for the purification of gangliosides from brain tissue. In Brown

F, ed. Transmissible Spongiform Encephalopathies - Impact on Animal and Human Health. Basel, Switzerland: Karger; 1993: 187-194.

Rosenberg R. et al. Precautions in handling tissues, fluids, and other contaminated materials from patients with documented or suspected Creutzfeldt-Jakob disease. *Ann Neurol.* 1986; 19(1); 75 - 77.

European Standard EN 12442-1-3. Animal tissues and their derivatives utilized in the manufacture of medical devices. Analysis and management of risk. European Committee for Standardization. February 2000.

Peri-Strips Dry is a registered trademark of Synovis Life Technologies, Inc. (US Patent #5,752,965; Additional Patents Pending).

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EEC-Dokument. Guidelines for minimizing the risk of transmitting agents causing spongiform encephalopathy via medicinal products. Biologicals. 1992: 20; 155 - 158.

DiMartino A. et al. Inactivation of the scrapie agent in a scaled-down procedure for the purification of gangliosides from brain tissue. Brown F, Hrsg., Transmissible Spongiform Encephalopathies - Impact on Animal and Human Health. Basel, Schweiz: Karger; 1993: 187-194.

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DiMartino A. et al. Inactivation of the scrapie agent in a scaled-down procedure for the purification of gangliosides from brain tissue. In Brown F, ed. *Transmissible Spongiform Encephalopathies - Impact on Animal and Human Health*. Basel, Switzerland: Karger; 1993: 187-194.

Rosenberg R. et al. Precautions in handling tissues, fluids, and other contaminated materials from patients with documented or suspected Creutzfeldt-Jakob disease. *Ann Neurol*. 1986; 19(1); 75 - 77.

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## REFERENTIES:

World Health Organization. Public health issues related to animal and human spongiform encephalopathies: memorandum from a WHO meeting. Bul WHO. 1992: 70(2); 183 - 190.

EEC Document. Guidelines for minimizing the risk of transmitting agents causing spongiform encephalopathy via medicinal products. Biologicals. 1992; 20; 155 - 158.

DiMartino A. et al. Inactivation of the scrapie agent in a scaled-down procedure for the purification of gangliosides from brain tissue. In Brown F, ed. Transmissible Spongiform Encephalopathies - Impact on Animal and Human Health. Basel, Switzerland: Karger; 1993: 187-194.

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EU-dokument. Guidelines for minimizing the risk of transmitting agents causing spongiform encephalopathy via medicinal products. Biologicals. 1992; 20; 155 - 158.

DiMartino A. et al. Inactivation of the scrapie agent in a scaled-down procedure for the purification of gangliosides from brain tissue. In Brown F, ed. Transmissible Spongiform Encephalopathies - Impact on Animal and Human Health. Basel, Svejs: Karger; 1993: 187-194.

Rosenberg R. et al. Precautions in handling tissues, fluids, and other contaminated materials from patients with documented or suspected Creutzfeldt-Jakob disease. Ann Neurol. 1986; 19(1); 75 - 77.

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EEC-dokument. Guidelines for minimizing the risk of transmitting agents causing spongiform encephalopathy via medicinal products. Biologicals. 1992; 20: 155–158. 20; 155 - 158.

DiMartino A et al. Inactivation of the scrapie agent in a scaled-down procedure for the purification of gangliosides from brain tissue. I Brown F,

utgåva Transmissible Spongiform Encephalopathies – Impact on Animal and Human Health. Basel, Schweiz: Karger; 1993: 187-194.

Rosenberg R et al. Precautions in handling tissues, fluids, and other contaminated materials from patients with documented or suspected Creutzfeldt-Jakob disease. Ann Neurol. 1986; 19(1); 75 - 77.

Europeisk standard EN 12442-1-3. Animal tissues and their derivatives utilized in the manufacture of medical devices. Analysis and management of risk. European Committee for Standardization. Februari 2000.

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*Authorized Representative in the European Economic Area:  
AR-MED Limited (Company Number 3154742)  
Runnymede Malthouse  
Egham, Surrey TW20 9BD  
United Kingdom*

Surgical Innovations

A Division of Synovis Life Technologies, Inc.

2575 University Ave. W.  
St. Paul, MN 55114-1024 USA  
651.796.7300 800.255.4018  
651.642.9018 (fax)  
[synovissurgical.com](http://synovissurgical.com)

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Peri-Strips Dry ett registrerat varumärke som ägs av Synovis Life Technologies, Inc.